

MINI SCREW'S INDICATIONS

- ★ **Correction of Borderline Orthognathic Cases into Non Surgical Way**
- ★ **Vertical Dimension**
 - Intrusion
 - Deep Bite Correction
 - Correction of Mutilated cases
Ex. Molar Intrusion
 - Molar Uprighting
 - Orthodontic Disimpaction
Ex. Canine / 3rd Molar
 - Correction of Occlusal Cant
- ★ **Antero-Posterior Dimension**
 - Full Arch / Segmental Arch Distalization or Mesialization
 - Molar Mesialization / Distalization
- ★ **Horizontal Dimenstion**
 - Cross Bite / Scissor Bite Correction



TUSK DENTA

Tusk Denta
35 Bloom Drive,
Brampton, Ontario L6P 3C3
Canada
(E) - tuskdenta@gmail.com
(Web) - tuskdenta.com

For Enquiry in India
+91-9913642443 / +91-8128056983

TUSK DENTA
MINI SCREW

**REDEFINING
THE
BIOMECHANICS**



**YOUR PRACTICE
OUR QUALITY**



TUSK DENTA



Enlarged Neck

More space to hang close spring or power chain

Smooth Mushroom Head

For comfort & retention of elastic chain



2-Way Round Holes

For lever arm to solve impacted tooth

Diameter of 2mm

Resistance to fracture



Self Drilling Design



improving the efficiency of insertion to avoid any broken tip




Double Neck Design

Easy hygiene control & extra attachment

PLACEMENT SITES




CONSIDERATIONS



SIZES OF
TUSK DENTA
MINI SCREW

Ø2.0 x L10mm
Ø2.0 x L12mm
Ø2.0 x L14mm

Biomechanics of retraction between Tusk Denta Screw Vs Alveolar Screw



CREATIVE
INNOVATIVE
PRACTICE

EXPERIENCE
THE
DIFFERENCE